



## Colon Hydromat standard: for hydro-colon therapy

Colon Hydromat standard is a reliable hydro-colon therapy system to be used on an every-day basis in doctors' surgeries. Temperature pre-selection and a water flow counter enable therapists to regulate the treatment process to the largest possible extent.



### Highlights of this system:

- Highest safety thanks to tried and tested quality and state-of-the-art technology
- Precise notification about treatment status on a display
- Temperature pre-selection and water flow counter allow for individual settings
- Compact construction and reliable design for trouble-free every-day use
- User-friendly control elements guarantee for an uncomplicated treatment

### Technical data

Name:	Colon Hydromat standard
Power supply:	100-240 V~50-60 Hz
Water pressure (cold and hot):	2-4 bar
Size (WxHxD):	525 x 400 x 160 mm
Weight:	21.7 kg
Classification according to Directive 93/42 EEC	II A

■ Ord. no. 15800 left  
 ■ Ord. no. 15801 right