



## Habamat Aquaclean: efficient protection from bacterial reflux

Habamat Aquaclean is the ideal addition to the hydro-colon product portfolio by HERRMANN Apparatebau but it can also be easily combined with devices by other manufacturers. It reliably prevents bacterial reflux to drinking water pipes. Its construction complies with the European directive EN 1717 on the protection against the pollution of potable water installations. Thanks to its compact size, Aquaclean easily fits under every standard examination couch.



## Highlights of this system:

- Reliable protection from bacterial reflux
- Certified quality: free discharge according to EN 1717
- The system easily fits under every standard examination couch
- Aquaclean can be combined with all hydro-colon treatment systems on the market

## Technical data

Name:	Habamat Aquaclean
Water inlet pressure:	1.5 - 4 bar
Temperature warm water outlet:	adjustable from 35°C to 85°C
Water Flow:	max. 14 l/min
Voltage:	230 V / 50 Hz
Size (WxHxD):	430 x 550 x 540 mm
Weight:	32 kg
	■ Ord. no. 15508