



Medozon compact: a reliable system for normobaric ozone therapy

Excellent operational comfort and high safety - these are characteristics of systems by HERRMANN Apparatebau. Medozon compact has been developed according to our strict quality guidelines and can be used for almost the entire range of modern ozone treatment. The integrated colour touch screen allows for an uncomplicated handling.



Highlights of this system:

- Can be used for all types of normobaric ozone therapy
- High flexibility thanks to the possibility to draw up ozone into a syringe
- Can be used for the insufflation of the intestine
- High compatibility thanks to Luer Lock connections
- Multilingual operational guidance (German, English, Spanish, Italian, Turkish)
- Colour touch screen for a more comfortable handling

Technical data

Devicename:	Medozon compact
Power supply:	230V~/50 Hz oder 115V~/50 Hz
Input:	240 VA
Concentration ozone:	0,5-80 µg/ml
Flow through:	1.0 l/min
Dimensions (WxHxD):	424 x 185 x 388 mm
Weight:	13.1 kg
Classification according to Directive 93/42 EEC	II a

Exclusive for the Chinese market!