



## Medozon compact: a reliable system for normobaric ozone therapy

Excellent operational comfort and high safety - these are characteristics of systems by HERRMANN Apparatebau. Medozon compact has been developed according to our strict quality guidelines and can be used for almost the entire range of modern ozone treatment. The integrated touch screen allows for an uncomplicated handling.



## Highlights of this system:

- Can be used for all types of normobaric ozone therapy
- High flexibility thanks to the possibility to draw up ozone into a syringe
- Can be used for the insufflation of the intestine
- High compatibility thanks to Luer Lock connections
- Multilingual operational guidance (German, English, Spanish, Italian, Turkish)

## Technical data

Name:	Medozon compact	
Power supply:	230V~/50 Hz oder 115V~/50 Hz	
Input:	150 VA	
Concentration ozone:	2-80 µg/ml	
Flow through:	1.2 l/min	
Size (WxHxD):	424 x 185 x 388 mm	
Weight:	13.1 kg	
Classification according to Directive 93/42 EEC	II a	<span style="color: green;">■</span> Ord. no. 16608 <span style="color: grey;">■</span> Ord. no. 16609